Breezing Pro™ 技术认证*

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
Figure 1: Correlations between Breezing Pro™ and the Douglas Bag Method for (A) VO₂ and (B) VCO₂, and (C) REE for 44 subjects


Figure 2: Bland-Altman plot of calculated REE from N = 44 unique subjects at resting conditions, differential percent error $\left(\frac{[REE \text{ (by Breezing Pro } ^{TM}) - REE \text{ (by Douglas Bag Method)}]}{REE \text{ (by Douglas Bag Method)}} \times 100 \right)$ between Breezing Pro TM and Douglas Bag Method

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
References

- [1] W. D. McArdle, F. I. Katch, and V. L. Katch, "Exercise Physiology: Energy, Nutrition, & Human Performance," Lippincott Williams & Wilkins, 2007.
- [2] M. M. Manore, N. L. Meyer, and J. Thompson, "Sport Nutrition for Health and Performance," *Human Kinetics (Ed.)*, vol. Second Edition, 2009.

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
Validation Study: Wearable mask vs. Medical Graphics (Ultima system)

Validation Study: Wearable mask vs. Medical Graphics (Ultima system)

Collaborations with Dr. Doina Kulick and Dr. Bruce Johnson, Mayo Clinic Arizona

Note: series connection of the Breezing Pro and Medical Graphics systems is important to warrantee same breathing conditions are tested.

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- 1.Total of subjects tested: 90+ with Gold standard method and Breath-by-Breath reference metabolic cart.
- 2.The wearable device for VO₂, VCO₂, and Resting Energy Expenditure (EE) demonstrated between 90-100% accuracy with 95% confidence interval.